



QUS (TAI & TSI)

R20 Quick Guide



1. Start QUS



Acquire liver parenchyma without vessel and reverberation artifact.

B-mode image is placed in the right lobe of the liver at the intercostal space while the patient maintained breath-holding.

① QUS



Tab the [QUS] on the touch screen to start.

② Dual Live

You can select the [Dual Live] depending on your preference to display single/Dual. (2D on the left, QUS on the right)





2. Activate TAI™



1	TAI	Activate the [TAI] function.
2	ROI Position	The ROI box is placed in a relatively homogenous region of the right lobe of the liver at least 2cm below the liver capsule, while avoiding areas with focal fat sparing or deposition and large vessels.
3	Set	Press the [Set] button on the control panel to confirm the result, the result value (dB/cm/MHz, R2) is displayed on the screen.
4	Store with Result	Store image including the TAI/TSI results on the left of the screen.
5	Store When Set & Unfreeze	Store Image when press [Set & Unfreeze] on the contextual button. * It is a function that sets the results and saves the image simultaneously with unfreezing.

3. The result of TAI™



1 Result

- TAI™ (Tissue attenuation imaging) provides quantitative measurement of attenuation coefficients of tissues.
- R² (Measurement quality)
 - Higher value (R² > 0.6) indicates reliable result.
 - R² < 0.6 (Unreliable result) The Results in red color is recommended to avoid.

TAI [dB/cm/MHz]
0.75 (R² 0.57)


2 Scroll Result

Scroll the knob and select among the result and delete with push.

3 Color Map

Choose Color Map Type.



4 Alpha Blending



5 Unit

Change Unit between dB/m and dB/cm/MHz.

4. Activate TSI™



1 TSI	Activate the [TSI] function.
2 ROI Position	The ROI box is placed in a relatively homogenous region of the right lobe of the liver at least 2cm below the liver capsule, while avoiding areas with focal fat sparing or deposition and large vessels.
3 Set	Press the [Set] button on the control panel to confirm the result, the result value is displayed on the screen.
4 Store with Result	Store image including the TAI/TSI results on the left of the screen.
5 Store When Set & Unfreeze	Store Image when press [Set & Unfreeze] on the contextual button. * It is a function that sets the results and saves the image simultaneously with unfreezing.

5. Result / Report



1 Report

Select the [Report] button on the touch screen and the result value is saved on the report. You can check the mean and median value of consecutive measurements.

2 Result

All results of the liver analysis may be reviewed including TAI, TSI, Fat Fraction.

3 Fat Quantification

Reference guide for Fat Quantification

4 Fat Fraction

Shows TAI/TSI results in correlation with MRI-PDFF.

Patient ID	TAI,TSI	Name		
Exam Date	2025-05-16			
Abdomen				
Liver Analysis				
	(dB/cm/MHz)			
	TAI	TSI		
Mean	0.59	72.71		
SD	0.02	1.99		
Median	0.60	72.11		
IQR	0.03	3.89		
IR/Med.(%)	5.53	5.39		
Fat Fraction	< 3.00%			
Fat Quantification (TAI)				
Grade 0	Grade 1	Grade 2	Grade 3	
dB/cm/MHz	< 0.72	0.72 ~ 0.83	0.83 ~ 0.86	0.86 ≤
dB/m	< 252	252 ~ 291	291 ~ 301	301 ≤
Fat Quantification (TSI)				
Grade 0	Grade 1	Grade 2	Grade 3	
-	< 95.6	95.6 ~ 98.4	98.4 ~ 98.9	98.9 ≤

QUS

6. Assign TAI™/TSI™ to User Key



★ Assign TAI/TSI to User Key

Assign TAI/TSI to User Key to access QUS without additional keystroke.

- a User Key(Setup : Utility → Setup)
- b Customize UI-U5 from the set up.
- c Select item and set to assign.

- 1 TAI
- 2 TSI

Enter TAI or TSI mode with assigned User Key and set to get each result.

Mode change can be done with contextual key in the control panel.





[User/Contextual Key in the control panel]

User/Contextual Key Guide

- ★ Assigned key is shown in the monitor
- ★ Contextual [TSI] button in the monitor
- ★ Contextual [TAI] button in the monitor

7. TAI™ / TSI™ Guide



TAI™ / TSI™ Guide

- NPO for at least 4 hours.
- Intercostal scan, acquire image with horizontal liver capsule.
 - ✳ *To maximize intercostal space, put patient in the supine or left decubitus(30 degrees oblique to the left) and raise right arm.*
- ROI in the middle.
- Steady and slow breathing is required to get a stable result.
- Avoid Vessels, shadowing and place the ROI.
- Place ROI 3cm below Skin or 1.5cm below liver capsule to avoid Reverberation Artifact caused by Glisson's Capsule.
- ROI placement should be between 3cm-9cm of depth.
- Use default size of the ROI(2*4).

- The features, options may not be commercially available in some countries.
- Sales and shipments are effective only after the approval by the regulatory affairs. Please contact your local sales representative for further details.
- This Quick guide does not include all of the details of instruction, for more detail, please refer to R20 User Manual.
- Do not distribute this document to customers unless relevant regulatory and legal affairs officers approve such distribution.
- This User Quick Guide is based on R20 V1.01.
- Disclaimer: Some Images in this content were obtained from other system.

SAMSUNG MEDISON CO., LTD.

© 2025 Samsung Medison All Rights Reserved.

Samsung Medison reserves the right to modify the design, packaging, specifications, and features shown herein, without prior notice or obligation.