

# **FAQ**

## **How to Automatically Display Day of Cycle from LMP**

## [FAQ – Day of Cycle by LMP]


- Question  
Can the system automatically calculate the menstrual cycle day using LMP?

- Answer

■ Path: Setup > System > Patient

Set **Calculate Day of Cycle by LMP** to **On**.

- ▶ When this setting is enabled, the system will automatically calculate and display the **Day of Cycle** based on the entered **LMP** (Last Menstrual Period) in the patient registration screen.



- It helps streamline OB/GYN workflows by providing immediate cycle day information without manual calculation.